



Why Perform Chemistry in Flow?





Why perform chemistry in flow?



Overview

- Faster reactions
- Safer reactions
- Faster reaction optimisation
- Reaction conditions not possible in batch
- Reactions are usually more selective
- Scale up is easier in flow than batch
- Easy integration of reaction analysis
- Reactions are easier to work-up in flow



Explanations



How does flow chemistry achieve faster reactions?

- It is much easier to pressurize flow reactors
- Higher pressures enable higher temperatures
- Higher temperatures result in faster reaction rates

- Example:


$$k = A \exp \left[- \left(\frac{E_a}{RT} \right)^\beta \right]$$

- Syrris products can pressurize to 20bar
- This means the boiling point of solvents increase by 100-150C

Solvent	1 bar	7 bar	17 bar
Dichloromethane	41°C	109°C	153°C
Methanol	65°C	138°C	185°C
Water	100°C	181°C	231°C

- **The Arrhenius rate law tells us that reactions are 2 x faster for every 10C rise**
- **Therefore 100C rise = 2x2x2x2x2x2x2x2x2x2 faster (over 1000x faster)**

Example – super-heating in flow



In batch:

- Incomplete reaction using batch chemistry despite 1 week of reflux
- Yield of the reaction was poor.

In flow:

- It was possible to superheat the THF up to 140°C
- Optimisation: reaction time, reaction temperature and equivalents of reagents were varied for each of the two reactions
- 100% conversion achieved for a variety of substrates at 140°C with 1h residence time.




Why does flow chemistry enable safer reactions?

- The quantity of reaction occurring at any one time is minimised
- The surface area:volume ratio is 1000s times higher in a flow reactor
- Example:
 - If a 10L batch reactor explodes, this has serious consequences
 - The same 10L can be passed through a 10ml flow reactor, ensuring that only 10ml is reacting at any time
 - For a fast reaction e.g. 1 min reaction this only takes an overnight run
 - **In this case the risk is 1/1000!!!**

Example – safer reaction in flow

- Triazole synthesis at Wyeth



- Hazard: azide compound. Use in medium and large scale is prohibited
- β -azidoethyl phenyl sulfide: bp =65°C, TSU showed exotherm at 155°C
- “Continuous-flow reactions have the potential to be much safer than batch reactions, as only a small amount of reactive and potentially hazardous material is heated or converted to product at any given time.”




How can flow perform reaction optimization faster?

- In a flow reactor it is extremely easy to vary:
 - Reaction time
 - Vary total flow rate
 - Reaction temperature
 - Low thermal mass
 - The ratio of reagents
 - Vary flow rate ratio
 - Concentration
 - Vary solvent stream
- One reaction is flushed out by the next (separated by a solvent) therefore only one reactor is used.
- **This means 50-100 reaction conditions can be investigated with just 15mins set-up time**

Example – quick optimisation in flow

- Optimisation of thiazole synthesis at Burnham Institute

Table 1. Optimization of thiazole **2a** synthesis



entry	time (min)	temp ($^{\circ}$ C)	H ₂ O (equiv)	ratio (3:5a:2a) ^a
1	2.5	50	0	0 : 8.8 : 1
2	5.0	50	0	0 : 8.6 : 1
3	2.5	100	0	0 : 1 : 1.7
4	5.0	100	0	0 : 1 : 2.9
5	2.5	150	0	0 : 1 : 8.7
6	5.0	150	0	0 : 1 : 6.8
7	5.0	150	1.0	0 : 1 : 11.7
8	5.0	150	5.0	0 : 1 : 104
9	5.0	150	10	1 : 0 : 3.5

- Varying residence time, temperature and water equivalent
- Optimum conditions identified in 9 experiments and 37.5 min total experiment time




How can flow chemistry achieve reaction conditions not possible in batch?

- Mixing happens by diffusion
 - This is much, much faster and more reliable than in batch
- Because the reactors are pre-heated/pre-cooled, the reaction can change temperature almost instantly
- Example 1
 - Heat up and cool down times are much faster than a microwave , therefore ultra hot, ultra fast reactions are easily possible
- Example 2
 - Deprotonate a substrate at low temperature, then add a nucleophile and instantly heat to a high temperature

How can flow reactions be more selective?

- Poor selectivity comes from variations in temperature, concentration and addition /stirring rates
- Due to a high surface area:volume ratio and diffusion mixing, flow chemistry offers:

- Excellent temperature control
- Minimal concentration gradient



- **Flow enables much better selectivity**

Example – cleaner reaction in flow

- Nitration of pyrazoles at AstraZeneca

Scheme 2. Temperature-dependent pyrazole nitration






Figure 4. Temperature population in a batch reactor vs a microreactor and impact on byproduct generation.²

- “In the case of a hazardous reaction, the risk of an incident is minimised as the accumulation of potentially dangerous intermediates is avoided.”


Example – better control in flow

- Reaction of CaCl_2 and Na_2CO_3 to synthesize CaCO_3 in batch and flow
- Exactly the same concentration, temperature and reaction/residence time
- Quality and reproducibility clearly much higher in flow than batch

Batch



Flow






Why is scale-up easier in flow chemistry?

- For 10x or 100x scale up
 - It is possible to flow for longer to make more product
 - E.g. from the same tap you can fill a cup or a bath
- For 1000x +
 - The fundamental principles of a higher surface area to volume ratio means that scaling up flow will reduce the heat transfer effect
 - The ability to use static mixers means that mixing is faster and more reproducible
- **Flow saves time and money scaling up reactions**

Scale Up Example



- Synthesis of Salicylaldehyde allyl ether
- Reaction used by Ley group to demonstrate their Micro Capillary Flow Disc Reactor and the FRX FLLEX (Flow Liquid-Liquid EXtraction)



- Published reaction conditions were room temperature in MeCN for 113mins
- Work up involved dilution with EtOAc and wash with $\text{HCl}_{(\text{aq})}$ using the FRX FLLEX

Optimisation and Scale Up

- Reaction was optimized and scaled to a 16ml reactor
 - Residence time reduced from 113min to 13.3min
 - Conversion remained high (97%)
- Reaction ran overnight



- **Result was 376g of Crude Product (at 97% yield from HPLC)**

How is reaction analysis made easier?


- To analyse multiple reactions in batch chemistry requires multiple probes (one per reactor).
- In flow chemistry, many different reactions can flow “under” the same probe.
 - The Sampler and Dilutor can automatically divert a tiny amount of each reaction for analysis



- In flow chemistry, the reaction automatically moves to the analytical system because it is flowing

Why is reaction work-up easier in flow?

- Batch chemistry relies on a separate operation to perform a work-up e.g. aqueous work-up, filtration or solid phase scavenging
- Because the reaction is already mobile flow chemistry offers in-line (integrated) work-up:
 - Flow liquid-liquid extraction (FLLEX)




- Solid phase reagents/scavengers/filtration



Liquid Liquid Extraction (FLLEX)


Step 1. Separation



- Diffusion of molecules (extraction) between slugs occurs rapidly, < 3 secs
- Separator, based on hydrophobic membrane, results in continuous flow of worked up product
- Separation of solvent pairs not typically thought of a practical e.g. THF and Aqueous

Syrris FLLEX as Copper Scavenging Unit

- Cu-catalyzed azide-alkyne cycloaddition reaction (CuAAC) (Click Chemistry)



Entry	Reactor type ^[b]	Cu catalyst (amount)	Yield [%] ^[c]
1	Cu tube	–	51
2	PFA tube	Cu(CH ₃ CN) ₄ BF ₄ (2.5 mol %)	91
3	PFA tube	[Cu(phen)(PPh ₃) ₂]NO ₃ (1 mol %)	92 (88) ^[d]
4	PFA tube	[Cu(phen)(PPh ₃) ₂]NO ₃ (0.5 mol %)	73

- CuAAC reaction is an ideal candidate for flow chemistry as high temperature required and intermediate azide formed pose an explosion hazard.





Table 2. Copper scavenging by extraction with aq. EDTA in copper(I)-catalyzed azide-alkyne cycloaddition.^[a]



Entry	Extraction ^[b]	Cu content ^[c] [ppm]
1	–	3156
2	flow (single extraction stage)	159 ± 9
3	1 × off-line	226 ± 51
4	2 × off-line	124 ± 4
5	3 × off-line	39 ± 1
6	flow (second extraction stage)	97 ± 3
7	flow (third extraction stage)	14 ± 1

[a] Reaction conditions: phenyl acetylene (0.225 M), phenyl azide (0.24 M), [Cu(phen)(PPh₃)₂]NO₃ (1 mol%), NMP, 180 °C. [b] For more details, see Supporting Information. [c] Measured by ICP–OES analysis of the product, average of two independent experiments.

RESULTS:

- The copper content could be effectively reduced to 119 ppm in a single equilibrium extraction stage using Reaction flow/aq EDTA/EtOAc = 1/7/5





Figure 5. Collection of the two phases after phase separation in continuous flow: aqueous phase containing EDTA, copper, and NMP (left vial), and the organic phase containing triazole (right vial).

Note: The amount of copper allowed in API is 15 ppm.

Example – work-up in flow

- 3-Hydroxymethylindoles synthesis by O'Shea et al.



Scheme 1. Multi-step strategy for the continuous-flow synthesis of 3-hydroxymethylindoles **5** and their conversion to **7** and **8** by acid-catalysed nucleophilic substitution.

- “Our current goal is to develop automated, sequentially performed homogeneous reactions with in-line continuous liquid–liquid extraction of the products”



Summary of flow chemistry benefits

- Faster reactions
- Safer reactions
- Faster reaction optimisation
- Reaction conditions not possible in batch
- Reactions are usually more selective
- Scale up is easier in flow than batch
- Easy integration of reaction analysis
- Reactions are easier to work-up in flow



Any questions?