## AH/AHS-N Range Semiautomatic Autoclaves

AH-21 N2

AHS-50 N

## AH/AHS-N Range

**N Class:** Suitable to sterilize: culture media, glass, liquids, plastic, metallic elements, waste bags, etc.

### Operating Principal

AH/AHS-N table-top range provides a single solution for the multiples sterilization needs on laboratories including liquids, culture media, biological waste, contaminated media, instruments, glassware and other laboratory items.

Load has to be placed into the vessel's baskets and, after manual filling with distilled water, device starts to heat up and purge until set temperature is reached. In case of solid objects, discharge can be manually forced during cooling stage to reduce cycle time.

### Digital Microprocessor

Digital microprocessor for an easy programming and parameters selection for sterilization process.

### Programs

Microprocessor has 10 programs, from P0 to P9, and the first four are pre-defined and protected.

| Program N° | Temperature (°C) | Pressure (bar) | Sterilization time (min.) | Sterilization mode |
|------------|------------------|----------------|---------------------------|--------------------|
| <b>P0</b>  | 115              | 0,6            | 60                        | SOLIDS-1           |
| <b>P1</b>  | 121              | 1,1            | 30                        | SOLIDS-1           |
| <b>P2</b>  | 133              | 2,1            | 20                        | SOLIDS-1           |
| <b>P3</b>  | 121              | 1,1            | 20                        | LIQUIDS-1          |

The rest of the programs are user-defined with the following parameters setting:

- Sterilization temperature
- Sterilization time
- Heart probe selection

## AH/AHS-N Range

### Sterilization Process

**Heating Stage** – Fast saturated steam production by means of powerful heating elements.

**Sterilization Control** – A PT-100 temperature probe controls the chamber temperature, with the possibility of installing another flexible PT-100 temperature probe for liquid programs.

**Cooling Stage** – In solid program, discharge can be manually forced during cooling stage to reduce the cycle time.

### Security

All safety standards are covered with safety valve, safety thermostat and safety pressure switch, door opening prevention system in case of overpressure, open door sensor and water level detector. EU Directive and Standards:

**2014/35/UE** *Low Voltage*

**2014/30/UE** *Electromagnetic Compatibility*

**2014/68/UE** *Pressure Equipment*

**EN-61010-1**

**EN-61010-2-040**

**EN-61326**

**AD 2000 Merkblatt**

### Communications and Documentation

**RS-232 Port** – For live monitoring of sterilization cycle directly to PC. It is needed the customized software SW7000 (optional).

## AH/AHS-N Range

**Integrated Printer** (optional) – For quality printouts after each cycle. Intervals can be adjustable. Printer has to be requested before ordering.


### Components supplied with the autoclave AH-21 N2

1. Tray support for four trays.
2. Three stainless steel perforated trays. Ref. BAH-21
3. Auxiliary tray.
4. A holding clamp.
5. 1.5 m flexible tube to drain condensation reservoir.

### Components supplied with the autoclaves AHS-50/75 N

1. Support for five trays.
2. Two perforated trays made of stainless steel. Ref. BAH-B
3. Auxiliary tray.
4. 1.5 m flexible tube to drain condensation reservoir.

## AH/AHS-N Range

### Technical Data

| Referencia<br>Reference<br>Référence | Volumen total / útil<br>Total volume / usable<br>Volume total / utile | Ø x H mm  | L x D x H mm     | Temp. °C | Presión máx. bar<br>Max. pressure bar<br>Pression max. bar | W    | Kg  | Rango tiempo<br>Time range<br>Rang de temps |
|--------------------------------------|---|-----------|------------------|----------|--|------|-----|---|
| <b>AH-21 N2</b>                      | 22/21   | 250 x 430 | 560 x 660 x 425  | 100-134  | 2,1  | 2000 | 55  | 0 - 500 h                                   |
| <b>AHS-50 N</b>                      | 55/50   | 400 x 400 | 805 x 805 x 650  | 100-134  | 2,1  | 2800 | 95  | 1' - 250'                                   |
| <b>AHS-75 N</b>                      | 79/75   | 400 x 600 | 805 x 1005 x 650 | 100-134  | 2,1  | 3200 | 110 | 1' - 250'                                   |

AHS-50/75-N models with closing wheel

Electrical connection for equipment of 21, 50 and 75 L: 230V (1P + N), 50/60 Hz.  
Possibility of modifying the voltage on request.

Audible and visible end-of-cycle sterilization indicator

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#### FLEXIBLE "HEART" PROBE PT-100

- Recommended for sterilizing liquids in flasks of 1 litre or bigger.
- The control is done from the microprocessor.
- This accessory should be installed in factory.

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#### THERMAL PAPER PRINTER\*

- It prints program number, cycle number, temperature, date and hour of the run, and error messages.
- This accessory should be installed in factory.

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Referencia: **ITS**

#### TABLE TOP PRINTER\*

- It prints program number, cycle number, temperature, date and hour of the run, and error messages.
- RS-232 connection.
- Not compatible with touch screen circuits.

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\* Serie AES, AE-DRY

#### PC CONNECTION SW7000


- Communication software between the equipment and the PC.
- Connection for RS-232.
- USB Adapter.
- Not compatible with touch screen circuits.

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\* Todos los modelos / All models / Tous les modèles

#### STERILIZATION CONTROL TAPE

- Chemical indicator for steam. The change of color indicates that the materials have been processed, without being a guarantee of proper sterilization. Class 1.
- To verify that sterilization is effective, additional methods are needed as biological indicators EN ISO 11138.
- Tape roll of 50 m x 19 mm

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