

## Riva MiniPress M II

### Single Punch Tablet Press

The RIVA Minipress MII is a Single Punch, off-center, Tablet Press. It is ideal for the production of small batches as well as for Pharmaceutical Research & Development. It is designed under strict GMP regulations.



Riva MiniPress M II

#### Easy Access

The smart design of the RIVA Minipress M II makes it easy to mount/dismount the Punch-Die set, operate the machine and clean it quickly and easily. On the other hand, the compression compartment is easily accessible and is closed by means of a detachable transparent acrylic panels equipped with safety sensors. Furthermore, the removable parts are equipped with quick release clamps mechanism for a quick disassembly for service and cleaning.



## Operating Principle

The gravity fill shoe performs an oscillating movement induced by the central crankshaft. It is connected to the hopper by a flexible tube. The pressure system is designed to stand up to 60 KN. The pressure is automatically adjusted according to the filling level of the die. This can be regulated by adjusting the upper punch penetration by means of the crankshaft. The filling volume of the die can be adjusted by regulating the weight cam at the lower central shaft of the punch holder. All contact parts are made of GMP Compliant materials and can be sterilized in accordance with GMP.

## Advantages

- » Small footprint
- » Good solution for R&D and small batch production
- » Easy to install & maintain
- » Corrosion resistant, medical & food grade contact parts

## Features

- » Compact
- » Smart design
- » GMP compliant

## Standard Scope of Supply

The RIVA MiniPress M II (part. no. 40-02000) comes with the following standard scope of supply:

- » Electronic Control Panel
- » Upper and Lower Punch Holder
- » Intermittent Oscillatory Gravity Feeder
- » Safety Sensor System in the Compression Compartment
- » Integrated Dust Nozzle

## Options

- » Instrumentation for pressure force indication
- » Documentation: IQ, OQ and user manual
- » Die-Punch Set of different sizes and shapes
- » Massive mobile stainless steel table
- » Tool Kit Box for installation, maintenance and service
- » Spare part kit

## Application Specifications

Parameter	Specification
<b>Operation principle</b>	Single Punch, Off-Center
<b>Process range</b>	<ul style="list-style-type: none"> <li>» Tableting of flowable granules/powders</li> <li>» R&amp;D, galenic development and low volume batch production</li> </ul>
<b>Feed characteristics</b>	Freely flowable granules and powders
<b>Range of applications</b>	Pharmaceutical, herbal & alternative medicine, agriculture, chemicals, detergents

## Tooling Specifications

Parameter	Specification
<b>Tooling</b>	“D” Special Riva Design
<b>Die Dimensions</b>	Ø = 38.10 mm; L = 23.81 mm
<b>Upper Punch</b>	Ø <sub>Max</sub> = 30.00 mm; L = 50.00 mm
<b>Lower Punch</b>	Ø <sub>Max</sub> = 30.00 mm; L = 67.50 mm
<b>Max. Tablet Diameter</b>	Round 24 mm; Shaped 10 x 25 mm
<b>Die Max. Load Depth</b>	18 mm
<b>Upper Punch Penetration</b>	3 – 14 mm

## Technical Specifications

Parameter	Specification
<b>Number of Stations</b>	One
<b>Safety Mechanism</b>	Panel sensors outside compression compartment
<b>Adjustment of Pressure and Tablet Height</b>	Automatic adjustment of pressure according to the filling level of the die
<b>Adjustment of Tablet Weight</b>	By regulating the weight cam via a hand wheel
<b>Speed Control</b>	Frequency inverter
<b>De-Duster, Vacuum Supply Nozzle</b>	Allows up to 45 m <sup>3</sup> /h
<b>Powder feeding</b>	Intermittent oscillatory gravity feeder
<b>Max. compression force</b>	60 KN (6 ton)
<b>Max. production rate</b>	6000 Tablet/h Depends on Feed Properties, Instrument and Process Parameters
<b>Power supply</b>	One single phase, 230 V
<b>Electric motor</b>	1.5 kW (2 HP)
<b>Instrument dimensions</b>	Approx. 830 x 720 x 400 mm (width x depth x height)
<b>Packaging dimensions</b>	Approx. 800 x 450 x 750 mm (width x depth x height)
<b>Net/ Gross weight (without sieve stacks)</b>	Approx. 280 / 435 kg

We reserve the right to make technical changes without any prior notice.