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Fuego SCS Art.-Nr. 8.000.000

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Windshield

made of stainless steel, for a
stable flame under conditions
of strong air flow

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Spray Protector
for more safety when handling
with pathogenic material

Art.-Nr./Art.-No.
6.000.200

The range:

Fuego SCS Art.-No. 8.000.000

with foot pedal, IR-Sensor and LC-Display
Adjustable IR-Sensor reaction distance, DoubleClick IR-Sensor (connectable)
4 standard-programs for IR-Sensor and foot pedal with count down display
SCS (Safety Control System) with BHC (Burner Head Control)
Removable and decomposable burner head
Tilt mechanism, right / left
Holding device for 3 inoculation loop holders
Nozzles for natural gas, propane / butane gas
Turbo flame
Wrench 17 mm for gas connection
Screwdriver for burner head and cover of burner shaft
Tubing connector with swivel nut
Power connection EURO (UK, USA, Japan)
Instruction manual and 2-year warranty

Technical Data:

Technology Microprocessor, LC- Display
Programs
IR-Sensor: Start-Stop with timer, 1 sec - 2 h
Auto-Off with timer, 1 sec - 2 h
Foot pedal: Standard (flame during pressed foot pedal)
Start-Stop with timer, 1 sec - 2 h

Safety features

Safety Control System (SCS)
with gas safety cut off: ignition and flame control
temperature monitor
burner head clogging
and assembly monitor (BHC)
automatic unit switch off, 1 - 90 min
residual heat display

Gas supply and consumption

Gas connection: 1/4" left with gas filter
Gas types: I12ELL3B/P: natural gas E/LL, 18 - 25 mbar
liquid gas, 20 - 50 mbar
Connected load: 100 g/h liquid gas
Continuous cartridge operation: CV 360 - 40 min, Express 444 - 50 min,
CG 1750 - 150 min, C 206 - 170 min,
CP 250 - 210 min, CV 470 - 370 min

Temperatures

Flame temperature: 1350 °C on liquid gas
1300 °C on natural gas (E)
Temperature threshold level: 1.2 kW liquid gas, 1kW natural gas

Electrical

Operating voltage: 10 V DC / 2 VA

Mechanical

Casing and operating controls: stainless steel / glass, UV and solvent resistant
Burner head: removable and decomposable, stainless steel
Cover of burner shaft: Ø 23 mm, with drains
Measurements (w x h x d): 103 x 49 x 130 mm
Weight: 700 g

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Safety Adapter

with integrated
pressure regulator for
gas cartridges Express 444

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Modern Design
without grooves and
rough edges

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Turboflamme
im Handumdrehen eine
extra harte und stabile
Flamme
Turbo Flame
an extremely firm and
consistent flame with only
a few flicks of the wrist

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Fuego^{SCS}

Safety Enhanced Laboratory Gas Burner

The Fuego SCS represents a new status among laboratory gas burners. With its innovative design and the extensive standard features this laboratory gas burner is suitable for all flame related applications in the laboratory.

Modern Design

The unique shape is eye-catching. The low profile body facilitates ergonomic operation; the streamlined design reduces air flow disturbances in a clean room workbench to a minimum. Moreover, the Fuego SCS is very easy to clean – a laboratory gas burner without grooves and rough edges. Clogging is not an issue. Spilled liquids simply escape through an integrated drain while solid substances can easily be removed by opening the burner shaft. It is even possible to dismantle the individual components of the removable burner head for in-depth cleaning.

Top Rate Safety

The Safety Control System SCS demonstrates state-of-the-art safety technology which constantly analyzes potential hazards and, if necessary, initiates safety measures, such as an interruption of the gas supply. In addition to the ignition and flame control function and the overheating protection, the new safety concept also features continuous burner head control (BHC). BHC unfailingly identifies burner head clogging by liquids or solid substances and ensures the correct assembly of the burner head.

Exceptional passive safety features: After prolonged use, the residual heat display will signal that the burner head is still hot in order to protect the user from burns. The adjustable automatic cut-off system inhibits unintentional ignition when the flame has not been ignited for a longer period of time.

Simple and Safe Operation

The Fuego SCS can be operated with a foot pedal or touch free IR-Sensor. For both options, different programs exist that have been developed from practical experience: flexible start-stop functions, continuous application for up to 2 hours or programs for short flame sterilization – exact to the second.

Moreover, the Fuego SCS automatically recognizes whether a foot pedal or the auto**loop**^{PRO} carousel for inoculation loops is connected, and displays other functions.

Convenience

The clear and easy to read display will always keep you in control. The logical multilingual menu facilitates the fast activation of any function. High-precision burning time regulation and IR-Sensor reaction distance adjustment are easily possible, among other options. When the burner is activated the display shows the remaining burning time.

Exclusive – The DoubleClick IR-Sensor

When activated, this additional safety function ensures that the burner can only be ignited by activating the DoubleClick IR-Sensor twice. As a result, unintentional ignition or ignition due to dropping or falling objects is virtually impossible, increasing safety considerably.

Another Plus

Excellent energy efficiency reduces the gas consumption. This increases the burning time when gas cartridges are used.

Robust

Fabricated entirely of stainless steel, with fireproof controls and a display, protected by heat resistant glass, the Fuego SCS can withstand extreme lab conditions.

Fuego SCS - Ideal for use in cleanroom workbenches and laboratories.