


VELP Scientifica Solutions for Overhead Stirring



ES

BS

DLS

LH

DLH

PW

Stirring speed rpm		from 50 to 1300	from 50 to 2000	from 50 to 2000	from 50 to 2000	from 50 to 2000	from 20 to 1200
Stirring volume L		up to 15	up to 25	up to 25	up to 40	up to 40	up to 70
Maximum viscosity mPa*s		1,000	10,000	25,000	50,000	50,000	100,000
Maximum torque Ncm		17	75	45	90	90	135
Chuck capacity mm		from 1 to 10	from 1 to 10	from 1 to 10	from 1 to 10	from 1 to 10	from 1 to 10
Counter reaction				•	•	•	
Dimensions (WxHxD) mm		80x160x200	80x215x196	80x215x196	80x230x196	80x230x196	80x230x196
	(in)	(3.1x6.3x7.9)	(3.1x8.5x7.7)	(3.1x8.5x7.7)	(3.1x9.0x7.7)	(3.1x9.0x7.7)	(3.1x9.0x7.7)
Weight	Kg	1.7	2.8	3.0	3.4	3.5	3.4
	(lb)	(3.7)	(6.2)	(6.6)	(7.5)	(7.7)	(7.5)
Power Supply		115 or 230 V	115 or 230 V	115 or 230 V	115 or 230 V	115 or 230 V	115 or 230 V
Power		35 W	60 W	60 W	120 W	120 W	120 W



www.velp.com